




01. CAE logo animation (2 sec)




02. Intro text about the product (CALL OUT) (3 sec)




03. Camera focuses in on the vaccine patch with needed accessories (3 sec)




04. Patient scans the QR code on the card (3 sec)




05. Patient sprays disinfectant (2 sec)




06. Patient rubs palms (2 sec)




07. Patient rubs/cleans the application area on the arm. (3 sec)




08. Patient opens applicator container. (3 sec)




09. Patient removes applicator (3 sec)




10. Patient removes the sticker from the applicator (3 sec)




11. Patient firmly places the applicator on previously disinfected area. (3 sec)




12. Patient presses the applicator at center wedging the sides upwards. (4 sec)




13. Patient waits for a few seconds as the vaccine is administered (TBD). (3 sec)




14. Patient removes the patch. (3 sec)




15. Patient disposes the patch back in it's container. (3 sec)




16. Patient applies pressure with a swab. (3 sec)



17. Patient sprays disinfectant once again. (3 sec)



18. Patient rubs hands once again. (3 sec)



19. Final product callout. (3 sec)



20. CAE logo animation. (2 sec)